

## BioNumerics Tutorial:

# Importing and processing VNTR capillary electrophoresis data

## 1 Preparing the database

1. Create a new database, install the *MLVA plugin* and create an MLVA scheme as described in the tutorial: "Setup of an MLVA scheme".


## 2 Importing curve files

A batch of Applied Biosystems curve files can be downloaded from the Applied Maths website (<http://www.applied-maths.com/download/sample-data>, click on "VNTR sequencer trace files"). These example files will be used to illustrate the import steps in this tutorial. Seven VNTR targets per sample were amplified: two pools were generated (**MP1** and **MP2**), containing four and three PCR products respectively, using the color dyes **6-FAM**, **VIC**, **NED**, and **PET**.

1. Select **File > Import...** (📁, **Ctrl+I**) to call the *Import* dialog box, choose **Import curves** under **Fingerprint type data** and press **<Import>**.
2. Browse to the downloaded and unzipped example data folder `VNTR sequencer trace files` and select all Applied Biosystems curves files (extension `.fsa`). Press **<Open>**.

The files are displayed in the *Input* wizard page and the default suggested **Fingerprint file name** is the folder name.

3. Change the name to **Batch1** and press **<Next>** (see Figure 1).




**Figure 1:** Selected curve files and fingerprint file name.

The way the information should be imported in the database can be specified with an import template. In the example data set, the dye name can be parsed from the file information and the sample and pool information can be parsed from the file names. A new import template needs to be defined:

4. Press the **<Create new>** button to call the *Import rules* dialog box (see Figure 2).

The *Import rules* dialog box lists the information present in the selected files as **Source**, their linked **Source type** and the **Destination** component they are associated with (currently all set to **<None>**).



**Figure 2:** Import rules.

5. Select **Curve dye** from the list, select **<Edit destination>** and select **Fingerprint dye** as corresponding field. Press **<OK>**.

Next, we will specify a new rule that links the part of the file name appearing before the underscore to the **Key** field.

6. Select **Name** from the list, select **<Edit destination>** and select **Key** as corresponding field. Press **<OK>**.
7. Check the option **Show advanced options**, make sure the last row is selected in the grid panel and press the **<Edit parsing>** button.
8. In the *Data parsing* dialog box, fill in following data parsing string: "[DATA]-\*". The asterisk will serve as wildcard.
9. Press the **<Preview>** button and press **<OK>** when the parsing is correct (see Figure 3).

Next, we will specify a new rule that links the part of the file name appearing after the underscore and before the hyphen ("-") to the **Pool** field.

10. Select **<Add rule>**, select **Name** under **File** as data source, press **<Next>**, and select **Fingerprint pool** as data destination. Press **<Next>** once more.
11. In the *Data parsing* dialog box, fill in following data parsing string: "\*-[DATA]-\*". The asterisk will serve as wildcard.




Figure 3: Parsing string.

12. Press the **<Preview>** button and press **<Next>** and **<Finish>** when the parsing is correct (see Figure 4).





Figure 4: Parsing string.

The *Import rules* dialog box should now look like Figure 2.

13. Press the **<Preview>** button to check the outcome of all defined import rules (see Figure 5).
14. Close the preview and press **<Next>** to go to the next step.

In the example data set, the LIZ channel contains the size standard (GeneScan 500 LIZ).


15. Make sure **LIZ** is selected as **Reference dye** and press **<Next>** and **<Finish>**.
16. Specify a template name, e.g. **Import VNTR curve files** and press **<OK>**.
17. Make sure **<Create new>** is selected as fingerprint type experiment, select the new template and press **<Next>**.



Nr.	Fingerprint pool	Fingerprint dye	Key
1	MP1	6-FAM	N2006-0003
2	MP1	VIC	N2006-0003
3	MP1	NED	N2006-0003
4	MP1	PET	N2006-0003
5	MP1	LIZ	N2006-0003
6	MP2	6-FAM	N2006-0003
7	MP2	VIC	N2006-0003
8	MP2	NED	N2006-0003
9	MP2	PET	N2006-0003
10	MP2	LIZ	N2006-0003
11	MP1	6-FAM	N2006-0004
12	MP1	VIC	N2006-0004
13	MP1	NED	N2006-0004
14	MP1	PET	N2006-0004
15	MP1	LIZ	N2006-0004
16	MP2	6-FAM	N2006-0004
17	MP2	VIC	N2006-0004
18	MP2	NED	N2006-0004
19	MP2	PET	N2006-0004

**Figure 5:** Preview of the import rules.

18. Specify a name for the new base fingerprint type experiment (e.g. **MLVA**) and press **<OK>**. Confirm the creation of the new experiment in the database.




Create new

Provide a name for the new fingerprint type:

OK Cancel

**Figure 6:** New base fingerprint type experiment.

A fingerprint type needs to be present in the database for each pool and dye combination. The names of these fingerprint types are composed of the base fingerprint type name, followed by the pool name, and the name of the dye. If one or more of these fingerprint types are not present in the database, a new dialog box pops up, listing the missing fingerprint types (see Figure 7).



New experiment types

These 10 new experiment types will be created:

- MLVAMP16-FAM
- MLVAMP1LIZ
- MLVAMP1NED
- MLVAMP1PET
- MLVAMP1VIC
- MLVAMP26-FAM

Do you want to continue?

Yes No

**Figure 7:** Missing fingerprint experiments.

19. Confirm the creation of the missing fingerprint type experiments.
20. Press **<Next>** to confirm the creation of new entries in the database.

21. Make sure *Open curve preprocessing window* is checked in the last step and press **<Finish>**.

For each dye checked in the *Dyes panel* of the *Import data* dialog box, a new fingerprint file is created, composed of the file name specified and the name of the dye (e.g. Batch1\_LIZ). These files are displayed in the *Fingerprint files* panel. The reference file is shown in the **Link** column. Double-clicking on a fingerprint file opens the *Fingerprint* window. If lane information was imported with the individual lanes, this information is displayed in the *Fingerprint information* panel.

The imported fingerprint lanes are linked to new entries in the database. The lanes are linked to the corresponding fingerprint "dye" type. The names of these fingerprint types are composed of the base fingerprint type name, followed by the pool name, and the name of the dye. The fingerprint type experiments are displayed in the *Experiment types* panel.

Entries for which fingerprint data was imported are selected in the database.

After data import, the *Main* window looks as in Figure 8.





Figure 8: The *Main* window after import of the data.


### 3 Processing curve files

When the option *Open curve preprocessing window* was checked in the last step of the import routine, the *Fingerprint curve processing* window opens when pressing the **<Finish>** button. All channels from the run are automatically loaded and displayed in the *Fingerprint curve processing* window.




The *Fingerprint curve processing* window can also be called from the *Main* window by highlighting one of the channels in the *Fingerprint files* panel and selecting **Open fingerprint data...** (🖋️). Alternatively, you can first open the *Fingerprint* window with **Edit > Open highlighted object...** (🖋️, Enter) and then select **File > Edit fingerprint data...** (🖋️).

1. Click on the  icon left of the data channels in the *Channels* panel.


The data channels are now hidden from the view and its icons are displayed as .

2. Use the zoom sliders on the left and on top to optimize the display of the fingerprint curves.

Since the raw chromatogram files have not undergone any preprocessing, normalization will have to be performed. This requires a *reference system* to be defined, based upon the marker peaks available in the reference dye.

3. Make sure the reference dye is the only dye visible in the upper panel (see Figure 9).
4. Select **Bands > Search reference bands...** (, **Ctrl+F**) to call the *Search reference bands* dialog box.
5. Leave the default settings unaltered and press **<OK>**.

The bands that fall within the specified criteria are marked with a solid line at the band's position (see Figure 9).




**Figure 9:** The *Fingerprint curve processing* window only displaying the reference dye.

6. To have a reference system automatically created based on a lane containing commercial size marker, first highlight a suitable lane and then select **References > Define size standard...**

This will display the *Size standard* dialog box, from which a size marker can be selected (see Figure 10).

In the example curve files, the GeneScan 500 LIZ size standard is used as reference, containing 16 bands with known molecular weight.

7. Select **GeneScan 500 LIZ** from the list, select **Pattern match** and press **<OK>** twice.
8. Save the data to the database with **File > Save** (, **Ctrl+S**).




Figure 10: Select marker from the list.

The software will automatically create the reference system and calibration curve for each of the fingerprint types. Since this allows the calculation of metrics information, a metrics scale now becomes available in the upper part of the *Fingerprints* panel.

Normalization is achieved by assigning bands in the reference channel to external reference positions.

9. To normalize a complete run at once, select **Normalization > Auto assign reference positions (all lanes)...** (🔍, **Ctrl+A**), leave all settings unaltered and press **<OK>**.
10. When the assignment of the marker bands to reference positions is made, the data can be shown in normalized mode with **Normalization > Show normalized view** (🔍, **Shift+N**).
11. Click on the 📁 icon left of the data channels in the *Channels* panel.
12. Click on the 📁 icon left of the LIZ reference channel.

The data channels are now shown and the reference channel is hidden from the view (see Figure 12).

13. Select **Bands > Search data bands...** (🔍, **Ctrl+Shift+F**) to call the *Search data bands* dialog box.
14. Check **Remove doublets**, **Remove shadow bands**, **Filter by fragment length** and specify a minimum length of 50. Press **<OK>**.
15. Save the changes and close the *Fingerprint curve processing* window.

When saving the data to the database with **File > Save** (💾, **Ctrl+S**) in the *Fingerprint curve processing* window, the software will automatically create the reference system and calibration curve for each of the fingerprint types. We can check this from the *Main* window:

16. Double-click on the base fingerprint type **MLVA** in the *Experiments* panel to open the *Fingerprint type* window.
17. In the *Fingerprint type* window, call **Settings > Edit reference system**, or double-click in the *R01* panel to call the *Fingerprint Reference system* window.

A calibration curve for the reference system **R01** of the base fingerprint type is displayed (see Figure 13).

18. Close the *Fingerprint Reference system* window and the *Fingerprint type* window.




Figure 11: Normalized view - reference dye.





Figure 12: Data channels.





**Figure 13:** The *Fingerprint Reference system* window, with a calibration curve displayed.